

Installing the Device Drivers for the BGA244

If you are having trouble installing the device drivers for the BGA244, go to the following link:

<https://ftdichip.com/drivers/>

There should be a link on the page labeled: **Windows Driver Installer (VCP & D2XX)**, please [click here](#).

That link downloads the zip file listed here: https://ftdichip.com/wp-content/uploads/2021/02/CDM21228_Setup.zip

This should have Windows D2XX and VCP (COM) drivers version 2.12.28.0, which is the latest version of the drivers they have available.