



- NOTES:
- 1) MATERIAL: ALUMINUM 1/4" STOCK
 - 2) BLACK ANODIZE SURFACE
 - 3) BREAK ALL EDGES 0.003
 - 4) TOL — ± 0.010

STANFORD RESEARCH SYSTEMS		BASE PLATE		
DRAWN	11-12-85	SIZE	FSCM NO.	DWG. NO. SR640-3
ISSUED	AWR	SCALE	2x	SRS# 7-00000-720 SHEET
				REV. D