



**Stanford Research Systems, Inc.**

1290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.

Phone: 408-744-9040 FAX: 408-744-9049

Email: [janied@thinkSRS.com](mailto:janied@thinkSRS.com) .. Website: [www.thinkSRS.com](http://www.thinkSRS.com)

---

***Volatility Statement for SRS BGA244 (Binary Gas Analyzer)***

The BGA244 contains two types of memory, as listed below:

**Non-Volatile Memory:**

256 Mbits Flash      Contains instrument's firmware, as well as user modifiable instrument settings.

To clear all user settings:

Recall "Factory" setting from the Store|Recall menu.

**Volatile Memory:**

256 Mbits SDRAM      Used for operation

To clear:

Cycle the power

Please feel free to contact SRS should you have any further questions.