

Certificate of Volatility				
Model: DC205	Part Number:	Manufacturer: Stanford Research Systems, Inc		
		Street Address: 1290-D Reamwood Ave		
		City: Sunnyvale	State:	Zip: 94089
Item Nomenclature:				
Volatile Memory				
Does the item contain volatile memory (i.e., memory whose contents are lost when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (SRAM, DRAM, etc.): SRAM	Size: 256kB	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: Program working memory	Process to Clear: remove power
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Non-Volatile Memory				
Does the item contain non-volatile memory (i.e., memory whose contents are retained when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (BBRAM, Flash, EEPROM, etc.): FLASH	Size: 1MB	User Modifiable: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Function: Firmware storage	Process to Clear: Programmed at factory only
Type (BBRAM, Flash, EEPROM, etc.): EEPROM	Size: 2MB	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: Stores calibration data and last configuration of instrument	Process to Clear: Two steps are required: (1) Press and hold the [Backspace] key while turning on power to clear the stored user configuration. (2) Press and hold the [Cancel] button while turning on power; then press [Enter/Start] to overwrite User Calibration data.
Type (BBRAM, Flash, EEPROM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Media				
Does the item contain media storage capability (i.e., removable or nonremovable disk drives, tape drives, memory cards, etc.)? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Additional Information:				
Vendor/Manufacturer Representative Information				
Name:	Title:	Office Phone:	Fax/Email:	