



## Stanford Research Systems, Inc.

1290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.

Phone: 408-744-9040 FAX: 408-744-9049

Email: [janied@thinkSRS.com](mailto:janied@thinkSRS.com) .. Website: [www.thinkSRS.com](http://www.thinkSRS.com)

---

### Volatility Statement for SRS DG535

The DG535 contains a Static Ram (3.3V lithium battery backup) which holds calibration bytes and instrument settings

Location	Type	Size	Function	User Modifiable?	Data Input Method	To clear User Settings
U403 On the top pc board	Static RAM	8KX8 100nS	Instrument Settings and Cal Bytes	Yes	Front Panel	Power on BSP * Or remove the battery (BT1)

\* With power off, hold backspace button and then power-on instrument (clear instrument settings)

Please feel free to contact SRS should you have any further questions.